Coolpack™

Prod. Nr. 134

f. Knie / Ellbogen klein

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Prod. Nr. 134

for knee / elbow small

The Coolpack Products contain a non-toxic gel, which can be cooled or heated many times. Cooling/heating can be relaxing and comforting and can lessen pain, as well as the symptoms of certain sports injuries, health problems or muscle tenseness.

COOLING: freezer or refrigerator

Your Coolpack may be used cold - for injuries, pain relief, cooling, etc. Put product in a plastic bag and cool down 1-2 hours in freezer. The gel insert will always stay flexible. Do not apply cold longer than 20 min.! Prolonged use only on the order of your physician! If you have an unstable circulatory system, consult a physician before cold use.

HEATING - MICROWAVE:

First Use: Heat your Coolpack in steps as follows (this will establish the standard heating time on your microwave - and will always stay the same). On the first use, note heating times:

- **Heat Coolpack for 30 sec.**
- Test Coolpack for desired heat (through thin cotton).
- If not warm enough, heat for additional 10 sec. until the desired temperature is reached.
- If the Coolpack starts to balloon **STOP MICROWAVE IMMEDIATELY.** This means the Coolpack has been overheated. Check for damage / shorten heating time on next use!

WARNING!! >>> Do not heat longer than a total of 30 sec. at 1100 W., 50 sec. at 800 W. or 70 sec. at 550 W. If overheated, gel insert and cover may be damaged. OVERHEATING WILL INVALIDATE OUR PRODUCT GUARANTEE OF 6 MONTH.

⇒ **!!! REHEATING: microwave - within 2-4 hours**

Heat the Coolpack **only** for 20-30 sec. Stop microwave if Coolpack balloons!!

Future Uses: Add up total heating time of first use and apply this time to heat your Coolpack.

Caution: When used with another microwave, proceed as in First Use.

HEATING: hot water

Boil a pot of water and take it off the heating plate. Then, take the gel inserts out from their cotton covers and place them for 5-10 minutes into hot water. Carefully take the gel insert from the water and replace it in the cover. Your Coolpack is now ready for use.

Please note: The Coolpack is not supposed to heat intensely, but is supposed to provide comfortable warmth. Do not overheat! Keep in mind, that products heated in a microwave oven only reach their final temperature a short time after the end of the cooking cycle.

For your information: Microwave ovens vary according to Model, Type and Power (Watt). With lower power the Coolpack will need a longer heating time, with higher power a shorter time. **Do not heat gel insert without it's cover.**

ATTENTION: - DO NOT REUSE A DAMAGED PRODUCT.

- DO NOT LET SMALL CHILDREN OR BABIES PLAY WITH A HOT OR COLD GEL PACK.
- DO NOT PLACE HOT OR COLD GEL PACK DIRECTLY ON SKIN

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