Face : Prod.Nr. 030

Face+: Prod.Nr. 033

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Biomark Therapeutic Face Masks are innovative and patented products for personal use, efficient and easy to use. They contain special non-toxic fills and can be warmed and cooled many times. Also for use in combination with Creams, Aromatherapy and Magnets.

MODELS: **Face:** covers the face from ear to ear and from chin to forehead.
Face+: as above, but also with prolongations to cover the throat.

INSTRUCTIONS:

COLD APPLICATIONS: Used to firm the skin, or after surgery, peeling, dermabrasion, etc. Keep the Biomark Face Masks in your refrigerator, for ready use. The face masks may also be cooled shortly in a freezer, but care needs to be taken to not use them too cool. Always place wet cotton face shield between skin and mask.

If your Biomark face mask is used with our M+ soft, unipolar magnets, place them under the cool mask (small ones next to your eyes against crowsfeet, larger ones over lower jowls). Magnets promote higher blood flow in the treated area, which influences availability of nutrients and plumps up local skin.

WARM APPLICATIONS: Used e.g. to open pores, to make creams and nutrients available to the skin, treatment with aromatherapy oils, or simple use of warmth with sinus problems, etc. Always place wetted cotton face shield between skin and mask.

Warm application: for aromatherapy or stuffy sinuses - warm under hot running water.

Higher temperature: the safest and easiest way is to bring water to a boil, then take it from the heating plate and place the face mask 3-5 minutes into the water, depending on temperature desired.

Warm in microwave: Biomark Face Masks may also be heated in the microwave - but carefully! Microwaves come with various powers. The heating times given below are variable and need to be adapted to your oven, but will then be the same for every use. **Always test temperature before you place the mask on your face (over wet cotton face shield).**

Instructions:

- Heat with highest power setting
- Heat 40 - 80 Sec. to reach desired temperature of 50-55°C =120-128°F (with smaller ovens fold the mask and heat each side 20-40 seconds)
- **Test for desired temperature**
- Clean face and throat, then place wet cotton face shield (smooth side to face) on your face and heated masks on top
- Relax and apply for 10 minutes

Adaptation of masks to the forehead: Cut within the seam to the stamped hole, then attach the included pieces of hook and tape to opposite ends of the cut. The closure may now be pulled and fitted to your forehead.

Cleaning: Reusable Biomark Therapeutic Face Masks may be cleaned with common cleansers, e.g. hot or warm water, soapy solutions, window cleaning solutions, alcohol, etc. The cotton face shields can be washed out by hand or in a washing machine.

Warranty: We guarantee Biomark Face Masks 6 months (for private use) against manufacturing faults. If an unexpected fault should occur, please send the faulty product with proof of purchase to one of the addresses below or contact your supplier:

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